

α^+ Solution™
Plant Genomic DNA Extraction Mini Kit

Kit Contents:

Cat. No:	HBPGK 004 (4 preps_sample)	HBPGK 050 (50 preps)	HBPGK 100 (100 preps)
PG1 Buffer	2.0 ml	25 ml	55 ml
PG2 Buffer	1.0 ml	8 ml	15 ml
PG3 Buffer (concentrate) ^a	1.5 ml	15 ml	30 ml
GW Buffer (concentrate) ^b	0.8 ml	13 ml	26 ml
Wash Buffer (concentrate) ^c	1.5 ml	15 ml	30 ml
Elution Buffer	1.5 ml	15 ml	30 ml
RNase A (lyophilized)	1.5 mg	22 mg	43 mg
Filter Column	4 pcs	50 pcs	100 pcs
PG Column	4 pcs	50 pcs	100 pcs
Collection Tube	8 pcs	100 pcs	200 pcs
User Manual	1	1	1
Preparation of PG1 Buffer, GW Buffer and Wash for first use:			
Ethanol for PG3 Buffer ^a	3 ml	30 ml	60 ml
Ethanol for GW Buffer ^b	1.0 ml	17 ml	34 ml
Ethanol for Wash Buffer ^c	6 ml	60 ml	120 ml









Specification:

Principle:	mini spin column (silica matrix)
Sample size:	1. Wet weight \leq 100 mg 2. Dry weight \leq 20 mg
Operation time:	30 ~ 60 minutes
Binding capacity:	up to 60 μ g total DNA / column
Expected yield:	5 ~ 40 μ g / prep
Column applicability:	centrifugation and vacuum
Minimum elution volume:	50 μ l

Important Notes:

1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
2. Check PG1 Buffer before use, Warm PG1 Buffer at 60°C for 5 min if any precipitate formed.
3. Preheat dry baths or water baths to 65°C before the operation.
4. **Add required ethanol (96-100%) to PG3 Buffer, GW and Wash before use.**
5. Store RNase A at -20°C upon receipt of kit. Add sterile ddH₂O to RNase A tube to make a 50 mg/ml stock solution. Vortex and make sure that RNase A has been completely dissolved. Store the stock solution at -20°C.
6. All centrifuge steps are done at full speed (~ 18,000 x g) in a microcentrifuge.

Brief procedure:

Sample Preparation		Grind 50 mg of wet weight (up to 100 mg) plant tissue or 20 mg dry weight of plant tissue with liquid nitrogen.
Lysis		<ol style="list-style-type: none"> 1. Add 400 μl of PG1 Buffer. 2. Add 8 μl of RNase A stock solution (50 mg/ml). 3. Incubate at room temperature for 2 minutes then at 65 °C for 10~20 minutes. 4. Add 130 μl of PG2 Buffer. 5. Incubate on ice for 5 min.
Filtration		<ol style="list-style-type: none"> 1. Transfer the entire mixture to the Filter Column. 2. Centrifuge at full speed (~ 18,000 x g) for 3 min. 3. Transfer the clarified lysate to a new microcentrifuge tube (not provided). 4. Add 1.5 volume of PG3 Buffer (ethanol added).
Binding		<ol style="list-style-type: none"> 1. Transfer up to 750 μl of the sample mixture to PG mini column. 2. Centrifuge at full speed (18,000 x g) for 1 min. 3. Repeat this step for the rest of the sample mixture.
Wash		<ol style="list-style-type: none"> 1. Add 400 μl of GW Buffer. Centrifuge at full speed (18,000 x g) for 30 seconds 2. Add 650 μl of Wash Buffer Centrifuge at full speed (18,000 x g) for 30 seconds. Repeat this step for one more washing.
Dry		Centrifuge at full speed (~18,000 x g) for an additional 3 min to dry the column.
Elution		Add 50~200 μl of preheated Elution Buffer. Centrifuge at full speed (18,000 x g) for 1 min.
Pure DNA		

General Protocol

Please Read Important Notes Before Starting Following Steps.

STEP	PROCEDURE
1 Sample preparation	Grind 50 mg of wet weight (up to 100 mg) plant tissue or 20 mg dry weight of plant tissue with liquid nitrogen to a fine powder and transfer to a new microcentrifuge tube (not provided).
2.1 Lysis	<ol style="list-style-type: none"> 1. Add 400 µl of PG1 Buffer and 8 µl of RNase A stock solution (50 mg/ml) to the tissue powder and vortex vigorously. 2. Incubate the mixture at room temperature for 2 minutes then at 65°C for 10~20 minutes and invert 2-3 times during incubation.
2.2 Lysis	Add 130 µl of PG2 Buffer to the mixture. Vortex to mix well and incubate the mixture on ice for 5 min.
3 Filtration	<ol style="list-style-type: none"> 1. Place a Filter Column to a Collection Tube and transfer the entire mixture from previous step to the Filter Column. Centrifuge the Filter Column at full speed (~ 18,000 x g) for 3 min. 2. Transfer the clarified lysate (supernatant) from the Collection Tube to a new microcentrifuge tube (not provided). 3. Discard used Filter Column and Collection Tube. And adjust the volume of clarified lysate. 4. Add 1.5 volume of PG3 Buffer (ethanol added) to the clarified lysate and mix well by pipetting.
4 DNA Binding	<ol style="list-style-type: none"> 1. Place a PG Column to a new Collection Tube and transfer up to 750 µl of the sample mixture from previous step carefully to the PG Column. 2. Centrifuge at full speed (18,000 x g) for 1 min. Discard the flow-through and place the PG Column back to the Collection Tube. 3. Repeat this step for the rest of the sample mixture.
5.1 Wash	<ol style="list-style-type: none"> 1. Add 400 µl of GW Buffer (ethanol added) to the PG Column. 2. Centrifuge at full speed (18,000 x g) for 30 seconds than discard the flow-through.

5.2 Wash	<ol style="list-style-type: none"> 1. Add 650 µl of Wash Buffer (ethanol added) to PG Column. Centrifuge at full speed (18,000 x g) for 30 seconds than discard the flow-through. 2. Repeat this step for one more washing.
6 Dry column	Centrifuge the PG Mini Column at full speed for an additional 3 min to dry the PG Mini Column.
7 DNA Elution	<ol style="list-style-type: none"> 1. Add 100 µl of preheated Elution Buffer or ddH₂O (pH 7.5-9.0) to the membrane of the PG Mini Column. Stand the PG Mini Column for 3 min. 2. Centrifuge at full speed for 2 min to elute DNA.