

PH METER MICROPROCESSOR HBPH-10



FEATURES :-

- pH Range -2.00 to 16.00 pH.
- pH Resolution +0.01pH.
- pH Accuracy +0.01 pH.
- mV Range 0 to +1999mV.
- mV Accuracy ± 1 mV
- Temperature Range 0 to 1000 C, Manual / Automatic
- Temperature Accuracy +10 C.
- Calibration Points 3 Points.
- pH Buffer Options NIST (pH 4.00/6.86/9.18) or (pH4.01/ 7.00/10.00)
- Temperature Compensation 0 to 1000 C, Manual or Automatic.
- Power Requirements DC 9V, using AC Adaptor 220 VAC, 50 Hz.
- Dimensions 210 (L) mmX205 (W) mm X 65 (H) mm.
- Weight 1.75kg. OF pH



sales@helixbiosciences.com



www.helixbiosciences.com



+91-11-41613622
1800-890-4669

REGISTERED OFFICE :-

A-7, SOUTH WEST DELHI, NARAINA INDUSTRIAL AREA PHASE 1, NEW DELHI - 110028

MANUFACTURING :-

PLOT NO. 360, IMT INDUSTRIAL ESTATE, SECTOR 68, FARIDABAD, HARYANA - 121004