

VORTEX MIXER



Introduction :-

Mini Mixer is used for vibration, blend and stirring of sample tissue, cell, bacteria liquid and chemical reagent in life science and physical and chemical analysis field.

Feature :-

Mini and compact appearance powerful function, sucker type under-chassis, super shockproof, ideal for high speed operation. Multiple mixing modes, touch-operation, continuous operation and timing set. Speed adjustable, and wide speed range, Stable low speed, powerful high speed. Shaking blocks convenient for replacement, stable and reliable, sturdy and durable eccentric bearing designed. Various tube holders are for optional (Suitable for Eppendorf Tubes).

Technical Parameters: :-

Model	MIX-25P Rotation Type	MIX-30S Digital display
Operation display pattern	Knob + Scale	Key + digital display
Orbit	3 mm	2 mm
Shaking mode	circle	circle
Operation mode	Continuous running/ Short Mix	Continuous running/ Short Mix



VORTEX MIXER



Technical Parameters: :-

Model	MIX-25P Rotation Type	MIX-30S Digital display
Speed Range	0 ~ 2500rpm	300 ~ 3000rpm
Timing Range	/	1s ~ 999min
Voltage	AC220V/AC110V, 50/60Hz	DC24V
Power	60W	30W
Optional Accessories	All	All
Fuse	250V,1A, Ø5x20	/

Accessories :-



sales@helixbiosciences.com



www.helixbiosciences.com



+91-11-41613622
1800-890-4669

REGISTERED OFFICE :-

A-7, SOUTH WEST DELHI, NARAINA INDUSTRIAL AREA PHASE 1, NEW DELHI - 110028

MANUFACTURING :-

PLOT NO. 360, IMT INDUSTRIAL ESTATE, SECTOR 68, FARIDABAD, HARYANA - 121004

VORTEX MIXER



Exchangeable blocks for Mix series :-

Cat#	Description
MX-A	Standard head. Suit for tubes and containers with dia.<30mm
MX-P	Standard platform to hold foam tube holder
MX-B	Convex pad. Use with MX-P. Suit for tubes and containers with dia.<99mm
MX-C	Foam tube holder. For 0.2ml tubes, 70 holes. Use with MX-P
MX-D	Foam tube holder. For 0.5ml tubes, 36 holes. Use with MX-P
MX-E	Foam tube holder. For 1.5/2.0ml tubes, 20 holes. Use with MX-P
MX-F	Foam tube holder. For 15ml tubes, 8 holes. Use with MX-P
MX-G	Foam tube holder. For 50ml tubes, 4 holes. Use with MX-P
MX-H	Tube holder platform. For 1.5ml/2.0ml tubes x 22
MX-I	Tube holder platform. For 15ml tubes x 4
MX-J	Tube holder platform. For 50ml tubes x 2
MX-K	Plate holder platform. For all standard plates
MX-L	Universal platform for tube stands MIX- A/B/C and all standard plates
MX-M	Customized



sales@helixbiosciences.com



www.helixbiosciences.com



+91-11-41613622
1800-890-4669

REGISTERED OFFICE :-

A-7, SOUTH WEST DELHI, NARAINA INDUSTRIAL AREA PHASE 1, NEW DELHI - 110028

MANUFACTURING :-

PLOT NO. 360, IMT INDUSTRIAL ESTATE, SECTOR 68, FARIDABAD, HARYANA - 121004