





# TAQ Cycler

-  **Faster Heating & Cooling Rate**
-  **More Durable Structure Design**
-  **Smart Phone App Available**
-  **User Friendly Interface**



## REMOTE OPERATION INTERFACE



Remote Operation Interface on Mobile Phone App

Operation Interface on Device



## PRODUCT PARAMETER

MODEL	TC100 - 96G
Sample Capacity	96 x 0.2 ml (PCR tube/Plate)
Temp. Range	4 ~ 105 °C (min. resolution: 0.1°C)
Max Heating Rate	6 °C/sec
Max Cooling Rate	4 °C/sec
Temp. Uniformity	≤±0.3 °C
Temp. Accuracy	≤ 0.2 °C
Hot-lid Temp. Range	30 °C ~ 110 °C
Gradient Temp. Difference Range	1 - 30 °C
Memory	≥ 2000 programs on board
Network Connections	USB / LAN / WIFI
Hot-lid Mode	One - time locking with handle
Mobile App Compatibility	iOS / Android
Power Supply	600W
Dimension	420mm x 269mm x 254mm (LxWxH)
Net Weight	13kg
<b>Offer Price</b>	<b>235,000.00</b>



**Intelligent Operation**






**Temperature Accuracy**



**Mobile App Sync**

**ENQUIRE NOW** 

 **1800 890 4669**  
 **011-41613622**  
 **011-45062522**